

FIGURE 1




FIGURE 2




FIGURE 3




FIGURE 4

Properties of Diversa Fluorescent Proteins

DVSACyan

Number of amino acids	227	253
Calculated subunit mass (kDa)	25.9	28.6
Total mass (kDa)	51.8	57.3
Excitation maximum (nm)	448(463)	487
Emission maximum (nm)	491	507
Quantum yield	0.76	0.61
Extinction coefficient ($M^{-1} \text{ cm}^{-1}$)	18,900	98,200

Figure 5

DIVERSA

DVSAGreen vs. Other GFPs

Excitation Maxima




Figure 6

DVSAGreen vs. Other GFPs




Figure 7

DVSA/Cyan vs. Other Blue/Cyan FPs

Excitation Maxima





Figure 8

Figure 9


DVSA Cyan vs. Other Blue/Cyan FPs

Emission Maxima



Excitation and Emission Spectra

Diversa Fluorescent Proteins



*Spectra normalized to the peak excitation and emission fluorescence for each protein

DIVERSA

Figure 10

DVSA Green protein is brighter than EGFP

Quantum yield Extinction coefficient
($M^{-1} \text{ cm}^{-1}$) Relative brightness*

wtGFP	0.77-0.80 ^{1,2}	21,600-27,600 ^{1,2}	1
EGFP	0.6-0.7 ^{3,4}	39,200-55,900 ^{3,4}	1.42-1.77
DVSA Green	0.79	90,000	2.7-3.6
AmCyan	0.24 ²	40,000 ²	0.43-0.58
DVSACyan	0.76	18,900	0.65-0.88

* Relative brightness (maximal extinction coefficient multiplied by quantum yield) as compared to wtGFP

¹ Taken from Heim and Tsien, Current Biology 1996



² Taken from Matz et al, Nature Biotechnology, 1999

³ Taken from Zimmer, Chemical Reviews, 2002

⁴ Taken from Remington, Nature Biotechnology, 2002

DIVERSA

Excitation and Emission Spectra



*Spectra normalized to the peak excitation and emission fluorescence for each protein

Figure 12

Relative Brightness

Comparison with Commercially Available Fluorescent Proteins			
	Quantum Yield	Extinction Coefficient (M ⁻¹ cm ⁻¹)	Relative Brightness*
Discovery Point ¹ (wtAvGFP)	0.61	98,200 ²	2.7-3.6
Wild type AvGFP	0.77-0.80 ^{1,2}	21,600-27,600 ^{1,2}	1
EGFP	0.6-0.7 ^{3,4}	39,200-55,900 ^{3,4}	1.42-1.77
pGlow	0.79 ³	30,000 ³	1.1-1.4
Discovery Point ¹ (AmCyan)	0.76	18,900 ²	0.65-0.88
AmCyan	0.24 ²	40,000 ²	0.43-0.58
ECFP	0.4 ³	32,500 ³	0.59-0.78
* Relative brightness: (maximal extinction coefficient multiplied by quantum yield) as compared to wtAvGFP, + Measured per chromophore			

1. Heim and Tsien. Current Biology 1996
2. Matz et al. Nature Biotechnology. 1999
3. Zimmer. Chemical Reviews. 2002
4. Remington. Nature Biotechnology. 2002

Figure 13

Summary

Summary of Diversa's DiscoveryPoint™ Fluorescent Proteins

	DiscoveryPoint™ Green-PP	DiscoveryPoint™ Cyan-PP
Excitation/Emission max (nm)	487/507	448(463)/491
Stoke's shift (nm)	20	43(28)
Maturation time	Within 1 hour	Within 1 hour
Quantum yield	0.61	0.76
Extinction coefficient (M ⁻¹ cm ⁻¹)	98,200	18,900
Thermostable to 80°C	Yes	Yes
# of amino acids	228	227
Calculated subunit mass (kDa)	26.0	25.9
Total mass (kDa) - dimers	52.0	51.8

Figure 14

DIVERSA

